

AVC Fume Hood Control System

Award-winning Intuitive Software Streamlines Installation



Manufactured in the USA.





The award-winning, Accutrol AVC Fume Hood Control System with Insight software provides an unprecedented evolution in fume hood control technology. It significantly simplifies the start-up, operation and modification for VAV fume hoods. The Accutrol AVC Fume Hood Control System combines the award-winning design of the ISO 9001:2015 certified, low pressure drop AccuValve® Airflow Control Valve with a "smart" fume hood display and a powerful and intuitive user interface. This innovative combination allows owners to easily set up and change the fume hood configuration without complicated manuals.

Award-Winning Building Automation



The Accutrol AVC Fume Hood Control System – Winner of the 2016 AHR Innovation Award for Building Automation


Innovative Design

The fume hood control system is comprised of the revolutionary AccuValve with an integral airflow measurement and control system, which includes a "smart" display mounted on the front of the fume hood. This unique measurement and control system, referred to as ePI®, provides pressure-independent performance using sophisticated smart electronics rather than complex mechanical components in the air stream.

Our free AVC Insight software connects to the AccuValve AVC airflow control valve through a standard USB connection located on the fume hood display. This provides technicians with quick access to easily change configuration and set points as well as set alarm points.

A technician can set up their fume hood in as little as 8 minutes by simply moving from screen to screen and entering the information when prompted. The AVC Insight software eliminates the need for the complicated manuals of the past.

The Accutrol AVC Fume Hood Control System





Easy Configuration & Verification of Sash Positions

Variable Air Volume (VAV) fume hood types consisting of vertical sash, horizontal sash or combination sashes are easily configured through the intuitive Accutrol Insight software. This eliminates the need for complicated menu driven instructions. The Insight graphical user interface provides users with a diagram of the specific type of fume hood they are configuring and prompts the user to enter the pertinent sash dimensions. (See Figure 2, below.)

Once the sash configuration is complete, Insight prompts the user to open and close the sashes to perform a sash calibration. To verify that the sash calibration has been successfully completed, Insight provides a fume hood diagram showing the real-time vertical and horizontal sash positions. (See Figure 3, below.)

Insight employs this same intuitive graphical approach for all configuration parameters.




Figure 2

Vertical, horizontal or combination sashes are easily configured through the intuitive Accutrol Insight software.

Applications

- University research and development
- Life science
- Pharmaceutical
- Biotech
- Healthcare laboratories
- Secondary school laboratories
- Retrofit applications

Features & Benefits

- Integral operation with low pressure drop AccuValve®
- Visual/Audible indication with Flow, Alarm and Purge
- Provides text alarm messages on display
- Purge set point is programmable
- Up to 4 user-defined setback sequences
- Optimized energy savings
- Simple installation and wiring
- Simplified set-up (no menus)
- Intuitive graphical user interface using our Insight software




Figure 3

To verify sash calibration, Insight provides a fume hood diagram showing the real-time vertical and horizontal sash positions.

Unmatched Value

The Accutrol AVC Fume Hood Control System is designed with the owner in mind. Unlike other complicated fume hood control systems that are difficult to set up, maintain and change, our system provides you with a simple, intuitive and sustainable product for your most critical airflow control applications.

The Accutrol AVC Fume Hood Control System provides:

- Precise, sustainable low energy airflow control
- Lowers cost of ownership from installation through lifetime operation
- Improved safety and ease of use for laboratory staff
- Native BACnet® for industry standard protocol

Specifications

FUME HOOD DISPLAY (FHM3 TOUCHSCREEN)

Display 3.2" touchscreen

Display Mounting Requires single-gang electrical box, min. depth: 1-7/8" (47.6mm)

Power Powered through connection to

Requirements AVC6000 Airflow Control Valve

Face Plate ABS Plastic Enclosure and

Mounting Plate

Dimensions 2.8"W x 4.6"H (71mm x 117mm)

Connectivity RJ45 cable connection to

AVC6000 Airflow Control Valve

Alarm Conditions Low Face Velocity, High Face

Velocity, Face Velocity, Deviation, Low Volume, High Volume,

Volume Deviation

(INSIGHT Software Configurable)

FUME HOOD DISPLAY (FHMI)

Display Back-lit, 2-line I" (25mm) LCD

Display Mounting Requires single-gang electrical box,

min. depth: I-7/8" (47.6mm)

Power Powered through connection Requirements to AVC Airflow Control Valve)

Face Plate Anodized aluminum with

antimicrobial pressure

sensitive overlay

Dimensions 2.8"W x 4.6"H

 $(71 \text{mm} \times 117 \text{mm})$

Connectivity RJ45 cable connection to AVC

Airflow Control Valve

Alarm Conditions Low Face Velocity, High Face
Velocity Face Velocity Deviation

Velocity, Face Velocity Deviation, Low Volume, High Volume, Volume Deviation (INSIGHT Software Configurable)







Additional Reference for the AVC Fume Hood Control System

AVC Fume Hood Control System Video Demystifying Fume Hood Control Systems

For a free live demonstration at your facility please contact us at demo@accutrolllc.com.

Your representative is:

